

Why use the Guardian Reader?

Although the tests may be read visually, the *Guardian Reader*TM will document your results and will alert you to positive tests that may be missed due to poor lighting or because of low sample concentrations. It will also prevent you from erroneously reporting positive results by “seeing” lines that are not there.

As you look at the examples below, remember that the intensity of the sample line (below each red arrow) is proportional to the concentration of the biological agent present in the sample. All tests will display a visible control line (below each black arrow) of the same intensity, to confirm that the test has been run correctly. If necessary, use Adobe Acrobat’s Zoom Feature to get a better look at the lines.

